



INSTALLATION QUALIFICATION PROCEDURES

PIPETMAX® 268





The enclosed Installation Qualification (IQ) procedures are aimed at end users who are implementing GLP-type requirements and can be readily incorporated into proprietary Standard Operating Procedures (SOPs).




1. Instrument Identification

Manufacturer

Gilson, Inc.
 3000 Parmenter Street
 PO Box 620027
 Middleton, WI 53562-0027 USA

| Instrument Name | Record Serial Number | Serial Number Locations (Compare to ensure that they match) |
|-----------------------------|--|---|
| PIPETMAX® 268 | | <ul style="list-style-type: none"> Shipping label on the shipping container Inside panel of the right support on the instrument |
| | Record Serial Number(s) Pipette heads are required, but ordered separately. Record the serial number(s) below, if applicable. | Serial Number Locations (Compare to ensure that they match) |
| MAX8x20 Pipette Head(s) | | <ul style="list-style-type: none"> Shipping envelope Document supplied with the pipette head Pipette head |
| MAX8x200 Pipette Head(s) | | <ul style="list-style-type: none"> Shipping envelope Document supplied with the pipette head Pipette head |
| MAX1000 Pipette Head(s) | | <ul style="list-style-type: none"> Shipping envelope Document supplied with the pipette head Pipette head |

| Instrument Name | Firmware Version | How to Obtain the Version |
|-----------------|--|--|
| PIPETMAX® 268 | XY system version: _____ Z system version: _____ Pipette system version: _____ | <ol style="list-style-type: none"> Use the supplied USB cable to make the connection between the USB device port on the PIPETMAX® 268 and a USB port on the supplied tablet or PC. Power the tablet on (TRILUTION® micro starts automatically on the tablet) or start TRILUTION® micro (click the Start button and then select All Programs > Gilson Applications > TRILUTION micro x.x > TRILUTION micro x.x.) Log in as an administrator. Select Home  to go to the Home screen. Select Settings. Select About. Select System info. <p>The version is displayed on the System information screen.</p> <p>If versions are not displayed, PIPETMAX and the software are not communicating; check your settings (including serial number) and cabling.</p> |



Supplier

Organization: _____

Address: _____

Phone Number: _____ FAX Number: _____

Name of Service Technician: _____

User

Organization: _____

Department: _____

Site (Room): _____

Date of Installation: _____

Primary Contact: _____

Name: _____ Phone: _____

E-Mail: _____



2. Pre-Installation

Instrument Description

PIPETMAX® 268 is a small footprint personal pipetting station used to perform repetitive, precise, pipetting actions.

For technical specifications, refer to the *PIPETMAX® 268 User's Guide*.

Unpacking

When unpacking the PIPETMAX, check the contents against the information listed under the **Unpacking** heading in **Chapter 1, Introduction**, of the *PIPETMAX® 268 User's Guide*.

Retain all packing material so the instrument may be shipped safely in the future.

Installation Site Requirements

Stable, level surface

Power Requirements

External Power Supply

Voltage Input

Frequency: 50 to 60 Hz

Voltage: 100–240V AC

Voltage Output

Voltage: 24V DC

Current rating: 6.25A, 150W

Temperature Range

5°C–40°C

Dimensions (W x D x H)

PIPETMAX® 268 with rotating cover

54.4 x 65.5 x 53.1 cm (21.4 x 25.8 x 20.9 in)

PIPETMAX® 268 with rotating cover installed on optional riser assembly for off-bed tip disposal

54.4 x 65.5 x 69.6 cm (21.4 x 25.8 x 27.4 in)

PIPETMAX® 268 without cover

50.8 x 64.3 x 49.5 cm (20 x 25.3 x 19.5 in)

PIPETMAX® 268 without cover installed on optional riser assembly for off-bed tip disposal

52.3 X 65 X 65.8 cm (20.6 x 25.6 x 25.9 in)

Documents

In addition to this document, the following are provided. For reference purposes, specify the storage location for each.

| | |
|-----------|---|
| LT3200USB | USB, PIPETMAX 268 REFERENCE DOCS _____ |
| LT255518 | SETUP GUIDE PIPETMAX 268 _____ |
| LT255521 | DECLARATION OF CONFORMITY PIPETMAX 268 _____ |
| LT310168 | QC CHECKLIST, PIPETMAX 268 _____ |
| LT310144 | ITEMS INCLUDED CHECKLIST, PIPETMAX 268 _____ |
| LT255522 | VALIDATION CERTIFICATE TRILUTION MICRO _____ |
| LT255535 | VALIDATION CERT TRILUTION MICRO v1.0 SP _____ |



3. Installation

Setup Guide

When setting up the PIPETMAX for the first time, follow the instructions provided in the *PIPETMAX® 268 Setup Guide*.

4. Acceptance

If no installation problems were detected or installation problems that were detected were resolved, have the local Gilson representative who installed the PIPETMAX (if applicable) provide the information requested below.

Gilson Representative:

Organization: _____

Signature: _____

Date: _____

End user was provided with a copy of this document.